












[www.remedimedical.com](http://www.remedimedical.com)





**CATALOGUE**




**2 0 2 0**

**LABORATORY**



**DISPOSABLES**





Remedi Laboratuvar Ürünleri Pazarlama Ltd. Şti. It is a company that continuously improves itself with its vast experience and practical business knowledge on plastic injection and plastic blow molding.

Remedi Laboratuvar Ürünleri Pazarlama Ltd. Şti. By meeting the needs of the main industries, it serves the sector by producing the required product from the design and molding to its production in a fast and high quality manner, as well as producing its own plastic laboratory consumables and medical consumables.


Believing in the understanding that “the most valuable element is human”, Remedi regularly organizes in-service trainings in order to increase the work experience of its personnel and to increase their corporate knowledge.

## Bacteriological Petri sample dishes

### Petri Sample dishes,

It's a flat round disposable plastic container that has a food in it and is used to produce creatures such as bacteria and fungi.



- Petri is produced from polystyrene raw materials.
- It contains a high level of transparency.
- It is suitable for aseptic or gamma sterile application.



#### Ø60x15mm Petri Sample dishes

Product code	Sterile	Barcode	Packing
001	R	8682009086005	18pcs/bag,50bag/parcel


Product code	Sterile	Barcode	Packing
002	A	8682009086012	18pcs/bag,50bag/parcel



#### Ø120x15mm Petri Sample dishes



Product code	Sterile	Barcode	Packing
011	R	8682009086081	20pcs/bag,12bag/parcel

Product code	Sterile	Barcode	Packing
012	A	8682009086098	20pcs/bag,12bag/parcel





#### Lockable Rodac Petri Sample dishes

Product code	Sterile	Barcode	Packing
004	R	8682009086036	24pcs/bag,45bag/parcel



#### Rodac Petri Sample dishes



Product code	Sterile	Barcode	Packing
003	R	8682009086029	24pcs/bag,45bag/parcel



### Ø90x14mm Petri Sample dishes

Product code	Sterile	Barcode	Packing
005	R	8682009086043	21pcs/bag,24bag/parcel

Product code	Sterile	Barcode	Packing
006	A	8682009086050	21pcs/bag,24bag/parcel



### Ø90x14mm Divided Petri Sample dishes

Product code	Sterile	Barcode	Packing
009	R	8682009086067	21pcs/bag,24bag/parcel

Product code	Sterile	Barcode	Packing
010	A	8682009086074	21pcs/bag,24bag/parcel

## 1ul and 10 ul Loops,


**It's a disposable product for collection, deploy and vaccination.**

- Smoth surface creates a sooth and soft line withaut damaging the gel surface.
- The design of the product makes it easy to do job.
- Flexibility or hardness may be given according to product raw materials.

STERILE - R




STERIL - R




### 10ul Loop

Product code	Barcode	Packing
014	8682009086111	10adet/poşet,1200poşet/koli

Sterile	Color	Raw material
R	Blue	PE




Body part



## 1ul Loop

Product code	Barcode	Packing
013	8682009086104	10adet/poşet,1200poşet/koli


Sterile	Color	Raw materials
R	Blue-Yellow	PE



## 16x110 TUBE,


**It's a type of tube that has a food in it used to produce creatures such as bacteria fungi.**

- It's diameter is 16mm height is 110 mm.
- It has a gasket and leak proof.
- It suitable for autoclave to be made at 121°



## 16X110 Tube

Product code	Raw material	Barcode	Packing
022	Polikarbonat	8682009086432	2000pcs/parcel



## 16mm Cover


Product code	Barcode	Gasket	Packing
024	8682009086456	YES	900pcs/bag,10pcs bag/ parcel





## Tissue Tracking Tape With Cover,

It is the preparation of tissues for microscopic examination and examination with tissue follow-up cassettes.

- It is produced from P.O.M raw material
- There is a writing area in the front area of 45 degrees.
- Its dimensions in the rectangular area are 5.40x 1.00 mm.
- There are 6 color variations.




## Coverless Tissue Tracking Tape



Color	Produce code	Barcode	Packing	
	White	053-1	8682009086357	500pcs/bag,10bag/parcel
	Yellow	053-2	8682009086357	500pcs/bag,10bag/parcel
	Green	053-3	8682009086357	500pcs/bag,10bag/parcel
	Blue	053-4	8682009086357	500pcs/bag,10bag/parcel
	Pink	053-5	8682009086357	500pcs/bag,10bag/parcel
	Orange	053-6	8682009086357	500pcs/bag,10bag/parcel

Color	Produce code	Barcode	Packing	
	White	052-1	8682009086364	250pcs/bag,10bag/parcel
	Yellow	052-2	8682009086364	250pcs/bag,10bag/parcel
	Green	052-3	8682009086364	250pcs/bag,10bag/parcel
	Blue	052-4	8682009086364	250pcs/bag,10bag/parcel
	Pink	052-5	8682009086364	250pcs/bag,10bag/parcel
	Orange	052-6	8682009086364	250pcs/bag,10bag/parcel



## Tissue Tracking Cassette With Cover



• Wide nail is designed for easy operation and easy to turn on and off.

• Lock system provides extra security during application.

## Compartmented tissue track tape




• It has a 45 write area.

• There is a separate biopsy area in itself

Color	Produce code	Barcode	Packing	
	White	051-1	8682009086371	250pcs/bag,10bag/parcel
	Yellow	051-2	8682009086371	250pcs/bag,10bag/parcel
	Green	051-3	8682009086371	250pcs/bag,10bag/parcel
	Blue	051-4	8682009086371	250pcs/bag,10bag/parcel
	Pink	051-5	8682009086371	250pcs/bag,10bag/parcel
	Orange	051-6	8682009086371	250pcs/bag,10bag/parcel

## 20 Lam Mape

Produce code	Raw materials	Barcode	Paking
049	Polystyrene	8682009086425	80pcs/parcel



## URINE CUPS

**It measure is between 20ml - 100ml.**

- Turn- open sealed special cover.
- It has a special area on which data can be written and label can be attached.
- It is produced with polypropylene raw materilas.

### 100ml Urine cups

Producte code	Sterile	Barcode	Packing
029	R	8682009086203	1pcs/1bag 150pcs/parcel

Producte code	Sterile	Barcode	Packing
028	A	8682009086197	1pcs/1bag 250pcs/parcel



### 100ml Urine cups

Producte code	Sterile	Barcode	Packing
027	Non-Sterile	8682009086180	750pcs/parcel



### 100ml Pneumatic Urine Container

**It is a lab material for use to urie container.**

- Turn- open sealed special cover.
- It has a special area on which data can be written and label can be attached.
- It is produced from polypropylene raw materilas..



### 100ml Pneumatic Urine Container

Producte code	Sterile	Barcode	Packing
035	A	8682009086418	1pcs/1bag 250pcs/parcel



### 100ml Pneumatic Urine Container

Producte code	Sterile	Barcode	Packing
034	Non-Steril	8682009086401	750pcs/parcel

## INTEGRATED URINE CONTAINER

**It measure is between 20ml - 100ml.**

- Turn- open sealed special cover.
- Suitable for vacuum system
- It is produced from polypropylene raw materials.



### 100ml Integrated Urine Container

Producte code	Sterile	Barcode	Packing
031	A	8682009086210	1pcs/1bag 250pcs/parcel



### 100ml Integrated Urine Container

Producte code	Sterile	Barcode	Packing
030	Non-Steril	8682009086227	300pcs/parcel

## 200ML URINE CUPS

**It measure is between 25ml - 200ml.**

- Turn- open sealed special cover.
- It has a special area on which data can be written and label can be attached.
- It is produced from polypropylene raw materials.

### 200ml Sample Container

Producte code	Sterile	Barcode	Packing
033	A	8682009086395	1pcs/1bag 150pcs/parcel



### 200ml Sample Container

Producte code	Sterile	Barcode	Packing
032	Non-Steril	8682009086388	300pcs/parcel



## WATER SAMPLE CONTAINERS

Water sample containers are a lab materials produced with thiosulfate and thiosulfate-free for samples taken from certain locations.

- It is divided into two with thiosulfate free thiosulfate.
- There are 2 different varieties opaque and amber colors
- The products are shipped as gamma sterile or non-sterile.



### 250ml Water Sample Container

Produce code	Sterile	Barcode	Tiosulfate	Packing
039	R	8682009086258	Yes	200pcs/ parcel

Produce code	Sterile	Barcode	Tiosulfate	Packing
037	R	8682009086234	-	200pcs/ parcel



### 250ml Water Sample Container

Produce code	Sterile	Barcode	Tiosulfate	Packing
040	R	8682009086265	Yes	200pcs/ parcel

Produce code	Sterile	Barcode	Tiosulfate	Packing
038	R	8682009086241	-	200pcs/ parcel



### 500ml Water Sample Container

Produce code	Sterile	Barcode	Tiosulfate	Packing
043	R	8682009086296	Yes	112pcs/ parcel

Produce code	Sterile	Barcode	Tiosulfate	Packing
041	R	8682009086272	-	112pcs/ parcel



## 500ml Water Sample Container

Produce code	Sterile	Barcode	Tiosulfate	Packing
044	R	8682009086302	Yes	112pcs/ parcel

Produce code	Sterile	Barcode	Tiosulfate	Packing
042	R	8682009086289	-	112pcs/ parcel



## 1000ml Water Sample Container

Produce code	Sterile	Barcode	Tiosulfate	Packing
047	R	8682009086333	Var	63pcs/ parcel

Produce code	Sterile	Barcode	Tiosulfate	Packing
045	R	8682009086319	-	63pcs/ parcel



## 1000ml Water Sample Container

Produce code	Sterile	Barcode	Tiosulfate	Packing
047	R	8682009086340	Yes	63pcs/ parcel

Produce code	Sterile	Barcode	Tiosulfate	Packing
045	R	8682009086326	-	63pcs/ parcel

## Cover with lock

Produce code	Raw materials	Cover color	Locked
YM-006	PE	White-opaque	Yes



- Size is Ø 45mm.
- Sealed.
- With foam gasket.



## 50ML FALCON TUBE

Used in clinical research labs such as daily cultural environment preparation and storage.

- It has a high transparent surface with flat cover.
- Easy to open- close cover.
- It is produced from PP material



- Special design applied for easy to use.
- It's indicative up 50ml.

- It has special are for write
- There are two types of skirts and conicals.



### 50ml Falcon tubes

Produce code	Cover color	Barcode	Packing
025	Blue	8682009086470	50/pcs,12/ parcel




### 50ml Falcon tube with skirt

Produce code	Cover color	Barcode	Packing
026	Blue	8682009086463	50/pcs ,12/parcel

### 2ml Collection tube

Produce code	Barcode	Raw materials	Packing
070	-	PP	100pcs/bag,100bag/ parcel





## 12x75 Tube

Produce code	Barcode	Raw materials	Packing
017	8682009086142	PS	500/bag,10/ parcel

Produce code	Barcode	Raw materials	Packing
018	8682009086159	PP	500/bag,10/ parcel



## 16x100 Round Tube

Produce code	Barcode	Raw materials	Packing
015	8682009086128	PS	250/bag,10/ parcel


Produce code	Barcode	Raw materials	Packing
016	8682009086135	PP	250/bag,10/ parcel




## 16x100 Conical Tube

Produce code	Barcode	Raw materials	Packing
019	8682009086166	PS	250/bag,10/ parcel

Produce code	Barcode	Raw materials	Packing
020	8682009086173	PP	250/bag,10/ parcel



• 16x100 Cover



• Easy to open- closed.

